





"Anatomy" of the GivMohr® Sling
 (see reverse for descriptions)



GivMohr Corporation
 3600 Osuna Road NE, Suite 316, Albuquerque, NM 87109
 Phone: 505-292-1144 Fax: 505-323-9526 Web: www.GivMohrSling.com





"Anatomy" of the GivMohr® Sling (see reverse for diagrams)

Handpiece

- Provides primary support of arm in reducing shoulder subluxation.
- OVAL - To maintain open web space. Can be modified to accommodate specific hand issues.
- CURVED - To support transverse arch of hand.
- Adjustments at handpiece:
 - To accommodate hand size.
 - To control for wrist flexion extension.

Elbow Strap

- Controls elbow position.
- Provides secondary arm support in reducing shoulder subluxation.
- Provides arm protection by keeping elbow close to body.

Anterior and Posterior Arm Straps

- Provide dynamic compressive forces through wrist, elbow and shoulder.
- Reduce potential for accommodation of nervous system.
- Protect shoulder from over-tightening (impingement).
- Allow arm swing with ambulation.
- Provide arm protection by keeping arm close to body.

Shoulder Harness

- Fashioned after "figure 8" harness.
- Facilitates thoracic extension, scapular retraction, and scapular depression.

Anterior and Posterior Buckles

- Adjust both buckles to reduce shoulder subluxation.
- Adjust posterior buckle independently once subluxation has been reduced to bring shoulder into neutral rotation.

Locking 4 Tab (back cross)

- Maintains proper positioning of shoulder harness.
- Position low on back for comfort. Lock in place.

Axillary Pad

- Position at anterior axillary wall of uninvolved arm for maximum comfort.

Wrist Cross

- Maintains positioning of sling to properly support arm.
- Helps hold wrist in neutral position.
- Position over dorsum of wrist.

