



GivMohr® SLING

Flaccid Upper Extremity Positioning Device
Designed & Developed by Victoria Givler OTR and Paul Mohr PT

Common Modifications to the GivMohr®

- [Adding wrist splint for proper wrist position](#)
- [Decreasing handle volume to accommodate distal motor return](#)
- [Mix sling sizes to accommodate different body types and arm lengths \(Bilateral\)](#)

Adding wrist splint for proper wrist position



Decreasing handle volume to accommodate distal motor return

Step #1

Cut plastic handle along webbing edge.

Step #2

Wrap handle with mole skin to reduce handle size and pad rough edges.

Mix sling sizes to accommodate different body types

Mixing shoulder harness and arm strap sizes is available to accommodate unique sizes and/or different arm lengths (bilateral sling) of individuals.

